

Useful tests for the differentiation of Group II *Enterococcus* species

Members of the genus *Enterococcus* are facultative anaerobic, catalase negative, gram positive cocci. The majority of the species grow in broth containing 6.5% NaCl, and they hydrolyze esculine in the presence of bile salts.

Group II *Enterococcus* can be differentiated from Groups I, III, IV, and V as they are acid formation from **Mannitol and Arginine positive**. The majority of species in this group have been recovered from human sources.

Species	Phenotypic characteristics					
	ARA	MOT	PIG	TEL	MGP	XYL
<i>E. faecium</i>	+	-	-	-	-	-
<i>E. casseliflavus</i>	+	+	+	-	+	+
<i>E. gallinarum</i>	+	+	-	-	+	+
<i>E. mundtii</i>	+	-	+	-	-	+
<i>E. lactis</i>	+	-	-	-	-	-
<i>E. faecalis</i>	-	-	-	+	-	-
<i>E. haemoperoxidus</i>	-	-	-	-	+	-
<i>E. rotai</i>	-	-	+	-	+	-
<i>E. thailandicus</i>	-	-	-	-	-	-

Names in bold: isolated from humans.

ARA: arabinose; MOT: motility; PIG: pigment; TEL: tellurite; MGP: methy- α -D-glucopyranside; XYL: xylose.

References:

Teixeira L.; Siqueira Carvalho M.; Fackam R.; Shewmaker P. *Enterococcus* in Manual of Clinical Microbiology. ed Jorgensen J.; Pfaller M.; Carrol K.; Funke G.; LandryM.; Richter S. Warnock D 2015, 11th edition. ASM Washington DC.

Identification of other *Streptococcus* Species: Streptococcus General Methods

<http://www.cdc.gov/streplab/strep-doc/index.html>