

VC172: VIDEO CHALLENGE QUESTIONNAIRE

Please follow this link and watch the video presented: <http://cmpt.ca/video-challenge-vc172/>
(Copy and paste the link to your URL).

After watching the video, answer the following questionnaire:

PART A - Indicate the source or error in the clip

- A. The suspicious colonies were not worked up.
- B. Chocolate agar was not included when setting up the sample culture.
- C. The swab submitted was not a rayon or nylon swab so it is not appropriate for this type of sample.
- D. Tellurite agar should be included when setting up this sample culture.

PART B – Select possible consequences of this error

- A. The culture results will be delayed at least 24 hours
- B. The pathogen may not be recovered if it is time/temperature sensitive or requires specific media.
- C. It may not be possible to collect a new sample before the patient is treated with antibiotics.
- D. The physician receives an incorrect report.
- E. A, B, and C are correct
- F. A, B, C, and D are correct

PART C – Check the proper action to remediate this error and reduce the probability of reoccurrence

- A. Ask Anne to complete an opportunity for improvement (OFI) report so that the laboratory error is recorded.
- B. Ensure that the sample set-up protocols are available and that staff are aware of and are trained for special set-up situations.
- C. Create a checklist to ensure set-up/accessioning staff includes appropriate media for specific clinical circumstances when setting up a sample.
- D. Ensure Anne receives extensive training on how to set up samples.
- E. A, B, and C are correct