CMPT Enteric Parasitology Quantitation Table

Quantitation		
Concentrate		
1+	Few	<11 per 22mm x 22mm coverslip
2+	Mod	11-20 per 22mm x 22mm coverslip
3+	Many	>20 per 22mm x 22mm coverslip
Stained smear		
1+	Few	<2 per 10 oil immersion fields (X1000)
2+	Mod	3-9 per 10 oil immersion fields (X1000)
3+	Many	>10 per 10 oil immerion fields (X1000)
Please note that quantitation is requested as an		
Internal check system for sample consistency		
and is not used as part of laboratory assessment		