

| W172 -July 10 , 2017: Membrane Filtration Method (cfu/100ml) - Time of Analysis (hours) |     |      |     |     |      |     |     |      |     |     |      |     |                     |                        |
|---|-----|------|-----|-----|------|-----|-----|------|-----|-----|------|-----|---------------------|------------------------|
| LabNo   | TC1 | ThC1 | EC1 | TC2 | ThC2 | EC2 | TC3 | ThC3 | EC3 | TC4 | ThC4 | EC4 | Time analysis (hrs) | TC media               |
| 1   | 42  | 47   | 42  | 17  | 0    | 0   | 28  | 28   | 28  | NT* | 0    | NT* | 24                  | mEndoles/mFC/NA<br>Mug |
| 1   | 42  | NT   | 42  | 18  | NT   | 0   | 19  | NT   | 19  | 16  | NT   | 0   | 24                  | CCA                    |
| 10  | 26  | NT   | 26  | 9   | NT   | 0   | 17  | NT   | 17  | 18  | NT   | 0   | 96                  |                        |
| 13  | 44  | 39   | 44  | 13  | 0    | 0   | 29  | 23   | 29  | 9   | 0    | 0   | 48                  |                        |
| 15  | 51  | 38   | NT  | 18  | 0    | NT  | 36  | 25   | NT  | 12  | 0    | NT  | 48                  |                        |
| 17  | 38  | NT   | 38  | 16  | NT   | 0   | 19  | NT   | 19  | 11  | NT   | 0   | 24                  |                        |
| 18  | 40  | 41   | 40  | 17  | 0    | 0   | 26  | 28   | 26  | 6   | 0    | 0   | 168                 |                        |
| 22  | 42  | NT   | 42  | 15  | NT   | 0   | 17  | NT   | 17  | 8   | NT   | 0   | 48                  |                        |
| 27  | 47  | 54   | 47  | 19  | 0    | 0   | 35  | 28   | 35  | 19  | 0    | 0   | 24                  |                        |
| 29  | 43  | 43   | 43  | 14  | 0    | 0   | 24  | 24   | 24  | 6   | 0    | 0   | 24                  |                        |
| 30  | 51  | 52   | 51  | 12  | 0    | 0   | 44  | 25   | 44  | 9   | 0    | 0   | 24                  |                        |
| 33  | 44  | 41   | 44  | 10  | 0    | 0   | 25  | 35   | 25  | 15  | 0    | 0   | 48                  | mEndoles/mFC/NA<br>Mug |
| 33  | 35  | NT   | 35  | 13  | NT   | 0   | 20  | NT   | 20  | 15  | NT   | 0   | 48                  | CCA                    |
| 36  | 49  | 48   | 49  | 12  | 0    | 0   | 31  | 22   | 31  | 10  | 0    | 0   | 48                  |                        |
| 47  | 49  | 47   | 49  | 12  | 0    | 0   | 23  | 25   | 25  | 15  | 0    | 0   | 24                  |                        |
| 48  | 37  | 46   | 37  | 15  | 0    | 0   | 27  | 26   | 27  | 11  | 0    | 0   | 72                  |                        |
| 61  | 45  | 51   | 45  | 19  | 0    | 0   | 21  | 19   | 21  | 13  | 0    | 0   | 24                  |                        |
| 80  | 48  | 45   | 48  | 16  | 0    | 0   | 27  | 26   | 27  | 16  | 0    | 0   | 24                  |                        |

\*- did not run test