

W183 - October 29, 2018: Membrane Filtration Method (cfu/100ml) - Time of Analysis (hours)

LabNo	TC1	ThC1	EC1	TC2	ThC2	EC2	TC3	ThC3	EC3	TC4	ThC4	EC4	Time analysis	Media
1	40	0	0	46	46	46	18	16	18	0	0	0	24	mEndoLes, mFC, Na MUG
1	35	NT	0	56	NT	56	13	NT	13	0	NT	0	24	Chromocult Agar
10	36	NT	<1	35	NT	35	18	NT	18	<1	NT	<1	72	
13	19	0	0	44	40	44	14	14	14	0	0	0	48	
15	25	0	NT	65	54	NT	23	13	NT	0	0	NT	48	
17	34	NT	0	45	NT	45	12	NT	12	0	NT	0	24	
18	39	0	0	40	62	40	16	15	16	0	0	0	72	
27	44	0	0	56	63	56	16	15	16	0	0	0	48	
29	30	0	0	45	45	45	13	13	13	0	0	0	24	
30	38	0	0	55	48	55	13	15	13	0	0	0	48	
33	36	0	0	46	42	46	18	16	18	0	0	0	360	mEndoLes, mFC, Na MUG
33	36	NT	0	52	NT	52	15	NT	15	0	NT	0	360	Chromocult Agar
36	28	0	0	46	56	56	13	17	17	0	0	0	24	
47	33	0	0	56	55	56	21	20	21	0	0	0	48	
48	35	0	0	60	49	60	13	22	13	0	0	0	264	
61	41	0	0	54	54	54	21	12	21	0	0	0	72	
80	34	<1	<1	67	67	67	21	21	21	<1	<1	<1	72	
95	7	NT	0	19	NT	19	15	NT	15	0	NT	0	144	