

**W202 - July 13, 2020: Membrane Filtration Method (cfu/100 ml) - Time of Analysis (hours)**

Lab no.	TC1	ThC1	EC1	TC2	ThC2	EC2	TC3	ThC3	EC3	TC4	ThC4	EC4	Analysis Time (hours)	Media
1	14	NT	0	0	NT	0	20	NT	0	33	NT	33	24	Chromocult
1	14	0	0	0	0	0	15	0	0	56	74	56	24	mEndoLes/mFC
10	16	NT	<1	<1	NT	<1	16	NT	<1	27	NT	27	96	
13	20	0	0	0	0	0	21	0	0	79	52	79	48	Chromocult
13	27	NT	0	0	NT	0	16	NT	0	62	NT	62	48	mEndoLes/NA Mug
15	18	0	NT	0	0	NT	17	0	NT	75	70	NT	48	
17	16	NT	0	0	NT	0	21	NT	0	50	NT	50	48	
22	15	NT	0	0	NT	0	13	NT	0	79	NT	79	72	mEndoLes/NA Mug
22	16	NT	0	0	NT	0	19	NT	0	56	NT	56	96	Chromocult
27w	12	0	0	0	0	0	17	0	0	53	70	53	48	
29	10	0	0	0	0	0	16	0	0	58	58	58	72	
30	20	0	0	0	0	0	14	0	0	65	65	65	336	
33	21	0	0	0	0	0	12	0	0	47	59	47	48	
36	25	0	0	0	0	0	12	0	0	48	46	46	72	
48	100*	<1	<2	<2	<1	<2	12	<1	<2	60	54	52	24	
61	14	0	0	0	0	0	17	0	0	75	54	75	48	
80	16	<1	<1	<1	<1	<1	17	<1	<1	72	72	72	24	
94	20	NT	0	0	NT	0	22	NT	0	51	NT	51	240	
95	1*	NT	1	0	NT	0	0*	NT	0	8*	NT	8*	360	

\*outliers not included in statistical calculations

NT - not tested